


RetCam Envision™ Connectivity Solutions

DICOM workflow option

RetCam Envision provides seamless integration and enhanced network capabilities. It's easier than ever to connect to the network, exchange and synchronize data.



Seamless connectivity with RetCam Envision's IT features include:

Encryption

Data can be encrypted on all RetCam Envision systems with Microsoft's Bitlocker disk encryption. RetCam Envision data files also support encryption, ensuring data security if transferred outside the RetCam Envision if transferred remotely.

Active Directory Integration

Use facility login credentials to secure access to the RetCam Envision system, eliminating local access control management.

Firewall & Antivirus

RetCam Envision systems are factory configured with firewall protection and are antivirus-enabled to maximize potential threats.

Audit Logging

Automated audit logs track user actions and system events to help manage security and compliance.

*NOTE: Not available in all countries. Please contact your local sales representative for availability.

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **033373 RevA**

natus.

Natus Medical Incorporated

natus.com